

PLASMA BASED STRUCTURED ILLUMINATION MICROSCOPY

*

In conventional wide-field fluorescent microscopy, a sample to be imaged is uniformly illuminated by light and the subsequent fluorescence is observed in the far-field through the objective of the microscope. The uniform nature of the illumination fundamentally restricts the resolution of the system to half the wavelength of light due to Abbe diffraction limit. With ever growing need to image tiny objects particularly in life sciences, modern microscopy techniques such as confocal and linear structured illumination microscopy use spatially non-uniform sources to illuminate the sample, resulting in resolution extending beyond the diffraction limit by a factor of 2 [?, ?]. In confocal microscopy, a small portion of sample is illuminated by a focused beam generated through a pinhole. The beam is laterally shifted to completely scan the sample, creating a sequence of images. Each image passes through another pinhole on the detector side. A high resolution image of the sample is generated by processing the image sequence, however confocal microscopy is a slow imaging technique. Moreover, part of light is discarded by the pinhole which may leave the signal strength from weakly fluorescent samples undetectably low. Structured Illumination microscopy (SIM) is a wide-field technique in which a fine illumination pattern such as a sinusoidal standing wave is used to generate *Moiré fringes* in the observed image. The high frequency content is mathematically reconstructed from a series of images acquired by shifting the pattern, yielding a high resolution image.

*Last Modified: 23:22, Monday 17th April, 2017.

1 Principle of Structured Illumination Microscopy

Consider $I(\mathbf{r})$ as the sinusoidal illumination intensity:

$$I(\mathbf{r}) = 1 + \cos(\mathbf{k}_\rho \cdot \mathbf{r} + \phi) \quad (1)$$

where $\mathbf{k}_\rho = k_x \hat{\mathbf{x}} + k_y \hat{\mathbf{y}}$ is the spatial frequency wavevector, $\mathbf{r} = x \hat{\mathbf{x}} + y \hat{\mathbf{y}}$ is the two-dimensional positional vector and ϕ is the pattern phase. The image of a sample $F(\mathbf{r})$ observed through a microscope can be expressed as:

$$M(\mathbf{r}) = [F(\mathbf{r}) \cdot I(\mathbf{r})] \otimes H(\mathbf{r}) \quad (2)$$

where $H(\mathbf{r})$ is the point spread function (PSF) of the microscope, and \cdot, \otimes denote multiplication and convolution operations in the spatial domain respectively. A frequency domain representation of the image by taking the Fourier transform is expressed as:

$$\begin{aligned} \tilde{M}(\mathbf{k}) &= [\tilde{F}(\mathbf{k}) \otimes \tilde{I}(\mathbf{k})] \cdot \tilde{H}(\mathbf{k}) \\ &= \frac{1}{2} [2\tilde{F}(\mathbf{k}) + \tilde{F}(\mathbf{k} - \mathbf{k}_\rho)e^{-j\phi} + \tilde{F}(\mathbf{k} + \mathbf{k}_\rho)e^{j\phi}] \cdot \tilde{H}(\mathbf{k}) \end{aligned} \quad (3)$$

where \sim over the letters indicates a frequency domain term and $\tilde{H}(k)$ is the optical transfer function (OTF) of the microscope. As evident in (3), a sinusoidal illumination pattern has three frequency components which generates an image which is linear combination of the sample along with two shifted versions as shown in Fig. 1(b). To reconstruct the sample, three different images need to be captured with different phase term ϕ . The process can be expressed as a system of linear equations,

$$\tilde{H}(\mathbf{k}) \cdot \begin{bmatrix} \tilde{F}(\mathbf{k}) \\ \tilde{F}(\mathbf{k} - \mathbf{k}_\rho) \\ \tilde{F}(\mathbf{k} + \mathbf{k}_\rho) \end{bmatrix} = \begin{bmatrix} 2 & e^{-j\phi_1} & e^{j\phi_1} \\ 2 & e^{-j\phi_2} & e^{j\phi_2} \\ 2 & e^{-j\phi_3} & e^{j\phi_3} \end{bmatrix}^{-1} \begin{bmatrix} \tilde{M}_1(\mathbf{k}) \\ \tilde{M}_2(\mathbf{k}) \\ \tilde{M}_3(\mathbf{k}) \end{bmatrix} \quad (4)$$

The phase shifts in (4) are known beforehand. Frequency content of the sample up to k_ρ can, therefore be observed due the Moiré effect which transports the high frequency information in to the observation region. To achieve two-dimensional enhancement in resolution, the sample has to be rotated about the optical axis of the microscope or the angular distribution of the illumination needs to be varied.




Figure 1: Resolution enhancement through SIM: (a) Diffraction limited observable region in frequency domain. Moiré effect using a sinusoidal illumination pattern bringing high frequency content under the observable region. The sample is rotated: (b) 0°, (c) 60°, (d) 120°. (e) Doubling of lateral resolution with effective coverage area twice the size of (a)

2 Generating Illumination pattern